Nitrile Gloves

DISPOSABLE GLOVES



MEDICAL GLOVES

Purpose

- To protect the patient and the professional
- To prevent microbes from getting attached to hands and spreading infection via hands
- To protect the hands from
- Occupational accidents
- Cleansing agents
- Chemicals
- Stick injuries

Use the gloves correctly

- Medical gloves are donned onto clean, dry hands
- · Hands are disinfected just prior to donning
- After disinfection hands must be allowed to dry before donning
- · Hands may get contaminated when undonning the gloves, therefore hands must be disinfected also immediately after undonning
- In all cases involving contaminating materials attention must be paid to correct undonning
- The outer surface of the glove must not get into contact with skin
- When necessary double donning can be applied -> decreased risk of skin contamination while undonning (e.g. when handling cytostatics)

Use of gloves

- Medical gloves are always single-use
- They are used per procedure and per patient
- Depending on the situation either sterile or factory clean gloves are used
- · Disposable gloves are neither washed nor disinfected



Nitrile gloves

Nitrile Examination Gloves

Comfortable, super soft flexible powder free nitrile gloves provide added safety in many applications. Designed with a special nitrile formulation they feel and fit like latex and allow full range of motion and excellent flexibility to minimize stress and fatigue. They do not contain natural rubber latex and are an excellent alternative for those suffering from Type I allergies.
















Features

- 100% latex free made from durable, puncture resistant, protein and powder-free nitrile, eliminates the Type I allergic reaction associated with natural rubber latex
- Highly elastic and super soft
- · Rolled rim to facilitate easy donning
- Designed to give a natural rubber-like feel
- Textured in finger tips for enhanced grip in all situations
- Food safe, suitable for handling fatty foods incl. cooking, olive and cod liver oil (Special Product)
- Protects against a wide variety of chemicals
- AQL 1.5, 2.5, 4.0 Available



- Complies with EN 455 and EN 374
- Complies with ASTM D6319 (USA Related Product)

Technical information

Туре	Powder-free, unsterile
Material	Nitrile butadiene rubber (NBR)
Colour	Blue, White, Purple, Pink, Black, Red
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm (≥230mm for USA)
Impermeability as per EN 455-1	AQL 1.5, 2.5, 4.0 Available
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years



Powder-free

REF	Size	Dispenser	Transport carton			
NG-141	XS	100	1000			
NG-141	S	100	1000			
NG-141	М	100	1000			
NG-141	L	100	1000			
NG-141	XL	100	1000			

Nitrile gloves

Nitrile Examination Gloves - Premium

On the market already for a long time nitrile gloves now have a new formula for increased softness and comfort and are even closer to latex gloves in strech, tactile sensitivity and performance. They are recommended for use in areas that may involve levels of stress.













Features

- 100% latex free made from durable nitrile, eliminates the Type I allergic reaction associated with natural rubber
- · Close to latex in stretch, tactile sensitivity and performance
- Designed to give a natural rubber-like feel
- Textured in finger tips for excellent grip in all situations
- Food safe, suitable for handling fatty foods incl. cooking, olive and cod liver oil. (Some Product)
- AQL 1.5

Quality standards

- Complies with EN 455 and EN 374
- Complies with ASTM D6319 (USA Related Product)
- Complies with ASTM F1671
- FDA 510(K) Available
- Approved to use with Chemotherapy Drugs

Technical information

Туре	Powder-free, unsterile
Material	Nitrile butadiene rubber (NBR)
Colour	Blue, White, Purple, Pink, Black, Red
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm (≥230mm for USA)
Impermeability as per EN 455-1	AQL 1.5
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years





Powder-free

REF	Size	Dispenser	Transport carton
NG-101	XS	100	1000
NG-101	S	100	1000
NG-101	М	100	1000
NG-101	L	100	1000
NG-101	XL	100	1000


Technical Data

Sizes	XS	S		М		L		XL	
Length (mm)	235±5mm	235±5mm		235±5m	m	235±5mm		235±5mm	
Palm Width	75±5mm	85±5mm		95±5mm	1	110±5mm		120±5mm	
Thumb L.	52±2mm	54±2mm		57±2mm	ı	59±2mm		65±2mm	
Forefinger L	66±2mm	68±2mm		72±2mm		74±2mm		80±2mm	
Middle finger L.	77±2mm	79±	±2mm	82±2mm		85±2mm		92±2mm	
Ring finger L.	68±2mm	71±2mm		75±2mm		77±2mm		82±2mm	
Little finger L.	52±2mm	54±2mm		56±2mm		60±2mm		64±2mm	
Dhariad		Viny l Premium		Synthetic Premium		Nitrile Premium		Latex Premium	
Physical Specification	Palm Thickness	Mir	n. 0.08mm	Min. 0.08mm		Min. 0.05mm		Min. 0.08mm	
	Finger Thickness		n. 0.05mm	Min. 0.05mm		Min. 0.05mm		Min. 0.08mm	
	Wrist Thickness	Mir	n. 0.06mm	Min. 0.06	imm	Min. 0.05mm		Min. 0.08mm	
	Tensile Strength Before Aging		n. 11MPa g. 12MPa	Min. 11M Avg. 16M		Min. 14MPa Avg. 18MPa		Min. 18MPa Avg. 20MPa	
	Tensile Strength After Aging		n. 11MPa g. 11MPa	Min. 11MPa Avg. 15MPa		Min. 14MPa Avg. 16MPa		Min. 14MPa Avg. 15MPa	
	Elongation Rate Before Aging	Min. 300% Avg. 420%		Min. 300% Avg. 460%		Min. 500% Avg. 580%		Min. 650% Avg. 680%	
	Elongation Rate After Aging		n. 300% Min. 300% g. 360% Avg. 410%			Min. 400% Avg. 470%		Min. 500% Avg. 600%	
Packaging	LxWxH		Method 1 (1000 PCs/CTN			d 2 PCs/CTN)		ethod 3 000 PCs/CTN)	
Dimensions	Dispenser		230x125x60mm		240x125x70mm		230x130x78mm		
(mm) & Loading Capabilities (Each product may differ from others)	Trans. CTN.	Trans. CTN.		315x258x250mm		385x260x255mm		408x268x240mm	
	Capability for 20 GP		1500 CTNs		1170 CTNs		1150 CTNs		
	Capability for 40 GP		3100 CTNs		2400 CTNs		2380 CTNs		
	Capability for 40 HQ		3450 CTNs		2700 CTNs		2600 CTNs		

The information above are for reference only. Specific order shall be confirmed according to regional regulations and customers' requiremts.

Qualification

MDD Product Qualification



FDA 510(K) EN 455 **ASTM** D5250 **ASTM** D6319

Quality Control System



GMP SGS EN ISO 13485 **SGS** EN ISO 13485 **SGS** EN ISO 9001

Other Qualification



日本食品衛生安全法検査成績書 Japanese Food Safety Code 370 BRC

CE-PPE-PVC

CE-PPE-NBR